Online Assessment (Even Sem/Part-I/Part-II Examinations 2019 - 2020

Course Name - Radiographic and Image processing Techniques Course Code - BMRIT204

- * You can submit the form ONLY ONCE.
- * Fill the following information for further process.
- * Required

1.	Email address *
2.	Name of the Student *
3.	Enter Full Student Code *
4.	Enter Roll No *
5.	Enter Registration No *
6.	Enter Course Code *

7. Enter Course Name *

8.	Select Your Programme *
	Mark only one oval.
	Diploma in Pharmacy
	Bachelor of Pharmacy
	B.TECH.(CSE)
	B.TECH.(ECE)
	BCA
	B.SC.(CS)
	B.SC.(BT)
	B.SC.(ANCS)
	B.SC.(HN)
	B.Sc.(MM)
	B.A.(MW)
	BBA
	B.COM
	B.A.(JMC)
	BBA(HM)
	BBA(LLB)
	B.OPTOMETRY
	B.SC.(MB)
	B.SC.(MLT)
	B.SC.(MRIT)
	B.SC.(PA)
	LLB
	PGDHM
	Dip.CSE
	Dip.ECE
	Dip.EE
	Dip.CE
	Dip.ME
	MCA
	() ()

	Offilite Assessment (Liver Seni/Fait-I/Fait-II Examinations 2019 - 2020
M.SC.(ANCS)	
M.SC.(MM)	
MBA	
M.SC.(BT)	
M.TECH(CSE)	
LLM	
M.A.(JMC)	
M.A.(ENG)	
M.SC.(MATH)	
M.SC.(MB)	

Answer all the questions. Each question carry one mark.

9. 1. Inherent film fog is

Mark only one oval.

- D = 0.2
- D = 0.12
- D = 0.4
- D = 0.5

10. 2. A screen converts

Mark only one oval.

- Energy of the light to X-ray beam
- Energy of the X-ray beam to light
- Energy of the light to UV rays
- Energy of the UV to light rays

11.	3. The substratum layer or binding layer is made-up
	Mark only one oval.
	Silver bromide crystals Gelatin plus acetone and water Gelatin None of these
12.	 4. The silver halide crystal used in emulsion of X-ray film is Mark only one oval. Silver bromide and Silver chloride Silver nitrate
	Silver iodide None of these
13.	5.Which of the following is not the advantage of double coated-ray film? Mark only one oval. With less exposure, better contrast Stability of film Prevent curling of film None of these

14.	6.Which of the following is not true double-coated film
	Mark only one oval.
	Crystal is 1 micrometer in size
	Each crystal can have 10 to the power 6 silver atoms
	Effected by mechanical pressure
	Emulsion is not sensitive to visible light
15.	7.The super coating is made-up of
	Mark only one oval.
	Cellulose
	Varnish
	Gelatin
	Silver iodide
16.	8.Single coated films are used for
	Mark only one oval.
	Mammography
	Ultrasound
	CT
	All of these

17.	9.The film boxes are arranged vertically upon the shelf in store room
	Mark only one oval.
	To avoid pressure on the film
	To avoid light leaking
	To minimize the space
	For convenience
18.	10.Unexposed silver halide has an electrostatic charge
	Mark only one oval.
	Positive
	Neutral
	Negative
	None
19.	11.Which is not true for green sensitive film?
	Mark only one oval.
	Less stable sensitometric response
	Improve image quality and sharpness
	Low replacement rates needed
	Better resistance to scratches

20.	12.vvnich is not true about film contrast?
	Mark only one oval.
	Inversely proportional to object film distance
	Inversely proportional to filtration
	Inversely proportional to film size
	Inversely proportional to focal spot size
01	42 The control of files decreased
21.	13.The speed of film depends on:
	Mark only one oval.
	Thickness of base
	Film Density
	Thickness of emulsion
	Latitude
22.	14.Scattered radiation can be minimized in PNS X-Ray by using of
	Mark only one oval.
	Cone
	Beam limiting device
	Filter
	Cylinder

23.	15.The primary function of the filtration is to reduce
	Mark only one oval.
	Image noise
	Scattered radiation
	Radiation dose
	To minimize the soft low energy photons-thus reducing patient skin dose.
24.	16.Which is not true about the optical density of the film –
	Mark only one oval.
	Density of the film is the sum of the densities on each side
	Inversely proportional to the number of the atoms in that area
	Directly proportional to the weight of the silver ions in that area
	Equal to the log 10 transmitted light- incident light
25.	17.The replenisher solution does not contain
	Mark only one oval.
	Metal
	Hydroquinone
	Phenidione
	Potassium bromide

26.	18.The working distance from the safe light to the film should not be less than
	Mark only one oval.
	0.10 meter
	1.2 meters
	5 meters
	15 meters
27.	19.Ceiling height in the dark room should be
	Mark only one oval.
	41 feet
	31 feet
	21 feet
	11 feet
28.	20.Amber safe light is compatible with films having
	Mark only one oval.
	Blue sensitive film
	Green sensitive film
	Yellow sensitive film
	None of these

29.	21. Which of the following material use for identification in radiography
	Mark only one oval.
	AI
	Cu
	CI
	Pb
30.	22.The first automatic film processor was developed in England without drying section by
	Mark only one oval.
	Kodak
	Siemens
	GE
	Potter
31.	23.The replenisher does not contain
	Mark only one oval.
	Potassium bromide
	Mettol
	Hydroquinone
	Phenidione

32.	24.The preservative in the developing solution is for
	Mark only one oval.
	To stop the development
	To decrease the rate of the development
	To prolong the effective life span of the developer
	None of these
33.	25.What does rinsing do?
	Mark only one oval.
	Remaining chemicals are washed from the film
	Removes excess developer from the emulsion; stops the development process
	Film is dried to remove the water and make it acceptable for handling and viewing
	Silver halide that was not exposed to radiation is dissolved and removed from the emulsion
34.	26.Single screen cassettes are used in
	Mark only one oval.
	Mammography
	Angiography
	Orthodontic radiography
	Ultrasound

35.	27.Intensifying screens are not used in
	Mark only one oval.
	Mammography
	Ultrasound
	Both Mammography and Ultrasound
	None of these
36.	28.The speed of the intensifying screen primarily depends upon
	Mark only one oval.
	Type of the film
	Size of the focal spot
	Size of the phosphor crystals
	Material used for the cassette
37.	29.Non-screen films are used for
	Mark only one oval.
	Thin part of the body
	Thick part of the body
	Less exposure technique
	None of these

38.	30.Fluorescence signifies
	Mark only one oval.
	Light emitted by screen within 10-8 sec
	Light emitted by screen after 10-8 sec
	Light emitted by screen at 10-8 sec
	None of these

This content is neither created nor endorsed by Google.

Google Forms